

## Powered stair climbing evacuation & rescue chair...

### What's an Escape-Chair® VOLT?

An Escape-Chair® VOLT is the **POWERED** solution to evacuate & rescue people to safety in a seated position up and/or down a flight of stairs. A simple procedure can make the Escape-Chair® VOLT fully functional within a few seconds.

### How does the Escape-Chair® VOLT work?

The Escape-Chair® VOLT is operated by a self-explanatory control panel with display on the operating handle. After unfolding the Escape-Chair® VOLT it is immediately ready for use. The single operator must push the operating handle down in the direction of the stairs.

### Advantages of the Escape-Chair® VOLT:

- Powered support makes going up the stairs much easier for operator.
- Ideal for quick evacuation (3 speed descent).
- Lithium cell enough for up/down 100 floors.
- Single operator chair.
- Ready for use within seconds.
- Easy to use and operate (wireless).
- Equipped with arm- and footrest.
- 3 height adjustable operating handle.
- Flame & heat resistant seat, headrest and dust cover.
- Battery indicator on display control panel.
- LED light for dark environment evacuation.
- CE certified.

### Our advice:


Make sure you have an Escape-Chair® VOLT available in your organization so that you're prepared for any evacuation and/or emergency, for up and/or down evacuation.

"Failure to prepare for an emergency has its consequences."





Display control panel



Powered tracks



Escape-Chair® VOLT folded



Heat and flame resistant  
Dust cover

| Specifications                                      | Escape-Chair® VOLT  |
|---|---------------------|
| Article number:                                     | 10050               |
| Dimensions folded up (height x width x depth):      | 115 x 50 x 27 cm    |
| Maximum load bearing:                               | 160 kg              |
| Maximum stair gradient:                             | 40°                 |
| Weight:   | ±30 kg              |
| Using time when full charged:                       | 2 hours             |
| Motor:  | 36V, 300W, no noise |
| Repeated charging times (on 220Volt):               | 500                 |
| Battery life:                                       | 3 - 5 years         |
| Extra batterypack:                                  | Optional            |
| Number of speeds:                                   | 3                   |
| Ready for use within seconds:                       | ✓                   |
| Foot rest:  | ✓                   |
| Arm rests:  | ✓                   |
| LED lightning for dark environment:                 | ✓                   |
| 3 height adjustable operating handle:               | ✓                   |
| Flame & heat resistant seat, headrest & dust cover: | ✓                   |
| Carrying handles front & back:                      | ✓                   |
| CE certified:                                       | ✓                   |
| Free demonstration and advise provided:             | ✓                   |
| Delivery incl. manual, identity sign & wall hooks:  | ✓                   |
| Warranty on non-wearing parts:                      | 1 year              |

Interim changes to technical specifications and printing errors reserved.



**Escape Mobility USA**

1261 S. Lyon Street, suite 405 • Santa Ana • CA 92705  
 T +1 888 258-0602 • F +1 714 912-2992  
 info@escape-mobilityusa.com  
 www.escape-mobilityusa.com